



## Key Locking

## DSC24SK-DP

### SHORT DILATOR HEPA DRYING CART for Savary-Gilliard® Dilators up to 28" in length

**Achieve the same drying capabilities of an endoscope drying cabinet with your Savary-Gilliard® Dilators!** This new medical innovation from the design and engineering team at Harloff allows you to vertically hang up to 24 dilators for channel drying, exterior drying and storage. This HEPA filtered dilator storage cart can hold dilators up to 28" in length, providing Positive Pressure savary dilator drying while removing 99.7% of particles, .03 Microns. The fan is adjustable from 10 to 180 minutes, and is set to activate each time the cabinet is accessed.



Engineered & Manufactured in the U.S.

Color shown is White with a Teal Door. More than 20 additional colors are available.

#### Features:

- Positive pressure fan system dries the full Savary-Gilliard® dilators and other wire guided dilators
- Includes 4 Savary pull out slides that each hold six dilators for easy access to similar sized dilators at one time
- Replaceable HEPA Filter removes 99.7% of particulates, .03 Microns
- Fan run time (from 10 to 180 minutes) is fully adjustable and activated each time the dilator storage cart is closed
- Durable scratch resistant heat sealed powder coat finish
- Each dilator storage cart can hold up to 24 individual dilators, varying in size from 5mm (15 French) to 20mm (60 French) and up to 28" in length
- Choice of lock options: latch lock, standard key lock (shown) or standard electronic lock
- Medical grade electrical system
- Sturdy, fully welded steel body
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Industry only fully curved protective bumper

#### Assembled Cart Dimensions:

- 40.38" H x 23.89" W x 22" D (102.6 cm H x 60.7 cm W x 55.9 cm D)
- Cart weight: 85 lbs (38.6 kg)
- Approximate shipping weight: 110 lbs (49.9 kg)
- Approximate Shipping Dims: 50" H x 29" W x 24" D (127 cm x 73.7 cm x 61 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/scope-drying-and-storage-cabinets/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300