



M-series

M3DS1830B05

M-SERIES TALL MINI WIDTH CART

with Five Drawers and a Breakaway Lock

Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 30" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and two 9" drawers (22.8 cm)
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Colors:

Shown in Red, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

All dimensions are available in the "MA-Dims-18wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Engineered &
Manufactured in
the U.S.



Tall Height, Mini Width, Five Drawer Emergency Crash Cart with a Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300