



M-series M3DS2424B05+MD24-EMG

M-SERIES SHORT EMERGENCY CART with MD24-EMG Emergency Package

Short, Medium Width, Five Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD24-EMG Accessories Include:



Oxygen tank
holder (O2HLDR)



IV pole (MD-IV2)



Cardiac board
and brackets
(MD-CARDBRD)

Color shown is
Red. More than
20 additional
colors are
available.



Engineered &
Manufactured
in the U.S.

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- 24" of vertical drawer space
- Cart is delivered assembled
- Removable plastic top
- Latex free
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Durable scratch resistant heat sealed powder coat finish
- Fully customizable drawer configurations



HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064
medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300

Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-24wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.