



**M-series**

**M3DS2430E06+MD24-ANS**

## M-SERIES TALL EMERGENCY CART with MD24-EMG Emergency Package

### Tall, Medium Width, Six Drawer Anesthesia Cart with Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### MD30-ANS Accessories Include:



3-Tier Back Rail  
(MD24-RBR3)



(2) 12" Side Rails  
(SIDERAIL12)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



Engineered &  
Manufactured  
in the U.S.

#### Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 30" of vertical drawer space
- Cart is delivered assembled
- Removable plastic top
- Latex free
- 3" (7.6 cm) single wheel stem caster set, full swivel - 2 with brake
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Durable scratch resistant heat sealed powder coat finish
- Fully customizable drawer configurations

Color shown is Navy.  
More than 20 additional  
colors are available.



**HARLOFF**  
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at  
<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064**  
**medicalsales@harloff.com**  
**www.harloff.com • Ph: (719) 637-0300**

#### Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-24wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.