



# M-series

## MDS1821E14

### M-SERIES X-SHORT MINI WIDTH CART

Four Drawers with an Electronic Pushbutton Lock

#### M-Series Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 21" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters
  - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

#### Colors:

Shown with a Beige body and Light Gray drawers, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).



Engineered & Manufactured in the U.S.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



### X-Short Height, Mini Width, Four Drawer Medical Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at <https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300