



M-series

MDS1824B05+MD18-PHB

M-SERIES SHORT PHLEBOTOMY CART

Mini Five Drawer Cart w/MD18-PHB Accessory Package

Short Height, Mini Width, Five Drawer Phlebotomy Cart with a Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Phlebotomy carts keep supplies nearby, are mobile, provide a dedicated work surface for procedures allow quick disposal of used needles. M-Series carts offer you choices.

MD18-PHB Accessories Include:



6-Compartment
Tilt Bin w/Rail Clips
(TILTBIN6RC)



Glove Box Holder
Single Direct Mount
GLOVEBOX1DM



Aluminum Waste
Container with Lid
AL680436



Push Handle
(MINI-PBHANDKIT)



3-Tier Back Rail
(MD18-RBR3)



Engineered &
Manufactured
in the U.S.

Color shown
in Light Blue.
More than
20 additional
colors are
available.



- Sturdy, uni-body steel construction
- Breakaway lock bar with a pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Fully customizable drawer configurations within 24" vertical space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters

HARLOFF

CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) & one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada.
1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document.
Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.