



# M-series

## MDS1824B05

### M-SERIES SHORT MINI WIDTH CART

Five Drawers with a Breakaway Lock

#### M-Series Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 24" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.

#### Colors:

Shown in Red, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



### Short Height, Mini Width, Five Drawer Emergency Crash Cart with a Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at <https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300