



M-series

MDS1824K03

M-SERIES SHORT MINI CART

Three Drawers with a Standard Key Lock

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 24" of vertical drawer space
- Drawer Configuration: One 6" drawer (15.2 cm) and two 9" drawers (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Colors:

Shown in Yellow, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Short Height, Mini Width, Three Drawer Infection Control Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



**Engineered &
Manufactured
in the U.S.**

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/isolation-infection-control/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**