



M-series

MDS1824K05+MD18-PHB

M-SERIES SHORT PHLEBOTOMY CART

Mini Five Drawer Cart w/MD18-PHB Accessory Package

Short Height, Mini Width, Five Drawer Phlebotomy Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Phlebotomy carts keep supplies nearby, are mobile, provide a dedicated work surface for procedures allow quick disposal of used needles. M-Series carts offer you choices.

MD18-PHB Accessories Include:



6-Compartment
Tilt Bin w/Rail Clips
(TILTBIN6RC)



Glove Box Holder
Single Direct Mount
GLOVEBOX1DM



Aluminum Waste
Container with Lid
AL680436



Push Handle
(MINI-PBHANDKIT)



3-Tier Back Rail
(MD18-RBR3)



Engineered &
Manufactured
in the U.S.

Color shown
in Hammertone
Blue. More than 20
additional colors
are available.



- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Fully customizable drawer configurations within 24" vertical space

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) & one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada.
1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.