



M-series

MDS1827B05

M-SERIES MEDIUM MINI WIDTH CART

with Five Drawers and a Breakaway Lock

Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 27" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters
 - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts


Colors:

Shown in Navy, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Engineered &
Manufactured
in the U.S.



Medium Height, Mini Width, Five Drawer Emergency Crash Cart with a Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300