



M-series

MDS1827WLP05

M-SERIES MINI WIDTH CART

with Wireless Electronic Lock & Prox Badge Reader

M-Series Features:

- Sturdy, uni-body steel construction
- Key lock override
- Membrane keypad with pass/fail and low battery light indicators
- Easily add new users with existing HID compatible badges (Prox and iClass)
- Wireless Lock can program up to 3,000 user codes and audit of the last 15,000 access attempts
- Enables you to use existing credentials. Allows for supervisor assignments; dual credential access available
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 27" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters
 - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Colors:

Shown in Mauve, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Medium Height, Mini Width, Five Drawer Medical Cart with a Wireless Electronic Lock and Prox Badge Reader

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

HARLOFF

CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300