



M-series

MDS1830B06+MD18-EMG

M-SERIES TALL EMERGENCY CART with MD18-EMG Emergency Package

Tall Mini Width, Six Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD18-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



Cardiac board and brackets (MD-CARDBRD)



Aluminum 16" Push Side Rail (SIDERAIL16)



IV pole (MD-IV2)



Engineered & Manufactured in the U.S.

Features:

- Sturdy, uni-body steel construction
- Durable scratch resistant heat sealed powder coat finish
- Breakaway lock bar or your choice of other lock options
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Removable plastic top
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters
 - one steering caster and two braking casters
- Body has pre-threaded accessory mounting holes
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Color shown in Red. More than 20 additional colors are available.

HARLOFF

CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

