



M-series

MDS1830E06+MD18-ANS

M-SERIES TALL ANESTHESIA CART

Mini Six Drawer Cart w/MD18-ANS Accessory Package

Tall Height, Mini Width, Six Drawer Anesthesia Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD18-ANS Accessories Include:



3-Tier Back Rail
(MD18-RBR3)



(3) 12" Side Rails
(SIDERAIL12)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



Engineered &
Manufactured
in the U.S.

Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Cart is delivered assembled
- Fully customizable drawer configs with 30" of vertical drawer space
- Industry best 5" quiet rolling casters - one steering caster and two braking casters

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Color shown
with a Light
Gray body and
Hammertone
Blue drawers.
More than 20
additional colors
are available.



HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.