



M-series

MDS1830L06

M-SERIES TALL MINI WIDTH CART

Six Drawers with a Basic Latch Lock

Features:

- Sturdy, uni-body steel construction
- Durable scratch resistant heat sealed powder coat finish
- Turn latch lock to secure drawers while moving cart
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 30" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Colors:

Shown in Hammettone Silver, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Engineered &
Manufactured
in the U.S.



Tall Height, Mini Width, Six Drawer Medical Cart with a Basic Latch Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

View similar products online at
<https://www.harloff.com/m-series-carts/>

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300