



MALIGNANT HYPERTHERMIA MDS2430B02-MH1

MALIGNANT HYPERTHERMIA CART with 1.0 Cubic Feet Medical Grade Refrigerator

Features:

- Low carbon powder-coated steel construction w/built-in Accucold® ARS1PV refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- High quality Accucold® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawers and refrigerator (refrigerator can be keylocked) with 100 breakaway seals
- Full extension drawer slides for easy access
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Cart ships ready to use with wheels installed
- Removable plastic top



Colors:

Shown in Navy, but available in more than twenty color options. View them online: www.harloff.com/finish-options/



Malignant Hyperthermia Cart, Two Drawer, Breakaway Lock with No Accessories

All dimensions are available in the "MA-Dims-24wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Engineered & Manufactured in the U.S.

Optional Accessories:

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Drawer Configuration:

- One 3" (7.6 cm) drawer, one 6" (15.2 cm) drawer & (1) Compartment with pre-installed refrigerator



Large wire basket



Side storage cabinets



Two large chart/document holders



Waste container & mounting bracket



View similar products online at <https://www.harloff.com/product-category/malignant-hyperthermia-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**