



# M-series

## MDS2430B06

### M-SERIES TALL EMERGENCY CART with Six Drawers and a Breakaway Lock

#### Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- 30" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawers (22.9 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-24wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

**Color shown is Red. More than 20 additional colors are available.**



### Tall, Medium Width, Six Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



**Engineered &  
Manufactured  
in the U.S.**

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300