



# MALIGNANT HYPERTHERMIA

**MDS2430K02-MH1**

## MALIGNANT HYPERTHERMIA CART

with 1.0 Cubic Feet Medical Grade Refrigerator

### Features:

- Low carbon powder-coated steel construction w/built-in Accucold® ARS1PV refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- High quality Accucold® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Standard key locks on drawers and refrigerator (refrigerator can be keylocked)
- Full extension drawer slides for easy access
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Cart ships ready to use with wheels installed
- Removable plastic top

**M-series**



### Colors:

Shown in Navy, but available in more than twenty color options. View them online: [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)



**Malignant Hyperthermia  
Cart, Two Drawer, Standard  
Key Lock with No Accessories**

**Engineered & Manufactured  
in the U.S.**

All dimensions are available in the "MA-Dims-24wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

#### Drawer Configuration:

- One 3" (7.6 cm) drawer, one 6" (15.2 cm) drawer & (1) Compartment with pre-installed refrigerator

### Optional Accessories:



Large wire basket



Side storage  
cabinets



Two large chart/  
document holders



Waste container &  
mounting bracket

**HARLOFF**  
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at  
<https://www.harloff.com/product-category/malignant-hyperthermia-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**