



M-series

MDS2430K06+MD24-ANS

M-SERIES TALL ANESTHESIA CART

Six Drawer Cart with MD24-ANS Accessory Package

Tall Height, Medium Width, Six Drawer Anesthesia Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD24-ANS Accessories Include:



3-Tier Back Rail
(MD24-RBR3)



(3) 12" Side Rails
(SIDERAIL12)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



Engineered &
Manufactured
in the U.S.

Features:

- Sturdy, uni-body steel construction with key lock and two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown with a Hammertone Blue body and Hammertone Gray drawers. More than 20 additional colors are available.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

All dimensions are available in the "MA-Dims-24wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.