



M-series MDS3024B05+MD30-EMG1

M-SERIES SHORT EMERGENCY CART with MD30-EMG1 Emergency Package

Short, Standard Width, 5-Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD30-EMG1 Accessories Include:



Drop shelf
(MD-DRPSHLF)



Oxygen tank
holder (O2HLDR)



IV pole
holder (MD-IV2)

Color shown is Red. More than
20 additional colors are available.



Cardiac board
and brackets
(MD-CARDBRD)



(2 Sets) Standard
3" Drawer Dividers
(MD30-DIV3-S)



(2) Utility hooks
(UHOOKS2)



Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Pack of 100 break-away seals
- 24" of vertical drawer space
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Removable plastic top
- Latex free
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters
 - one steering caster and two braking casters
- Heat sealed powder coat steel break-away lock bar
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish



Engineered &
Manufactured
in the U.S.

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

HARLOFF

CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com