



**M-series**

**MDS3024E14+MD30-ANS**

## M-SERIES SHORT ANESTHESIA CART

Four Drawer Cart with MD30-ANS Accessory Package

### Short Height, Standard Width, Four Drawer Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives.

### MD30-ANS Accessories Include:



3-Tier Back Rail  
(MD30-RBR3)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



(3) 12" Side Rails  
(SIDERAIL12)



Engineered &  
Manufactured  
in the U.S.

### Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles mounted on either side
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free

Color shown in Light Blue.  
More than 20 additional  
colors are available.



**HARLOFF**  
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

#### Drawer Configuration:

- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.5 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.