



M-series MDS3024K01-21DR

M-SERIES STANDARD WIDTH CART with One Drawer, Storage Area and Key Lock

Includes a 21" double door storage compartment (keyed separately) and adjustable interior shelf with six possible height locations

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top, easily replaceable
- 24" of vertical drawer space
- Drawer Configuration: One 3" drawer (7.6 cm) and one 21" double door storage compartment (53.3 cm)
- Latex free cart is delivered assembled
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Fully customizable drawer configurations
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Color shown is Navy. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Short, Standard Width, One Drawer Cart and Storage Compartment with Key Lock

Designed and built to meet the needs of today's healthcare market. With a side handle, comfort drawer pulls, expanded lock options and a standard plastic top to make things easier for your staff and to help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at <https://www.harloff.com/product-category/treatment-and-procedure-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**