



M-series

MDS3024K14+MD30-ANS

M-SERIES SHORT ANESTHESIA CART

Four Drawer Cart with MD30-ANS Accessory Package

Short Height, Standard Width, Four Drawer Anesthesia Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD30-ANS Accessories Include:



3-Tier Back Rail
(MD30-RBR3)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



(3) 12" Side Rails
(SIDERAIL12)

Color shown with a White body and Cream drawers. More than 20 additional colors are available.

Features:

- Sturdy, uni-body steel construction, delivered assembled
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles, mounted on either side
- 24" of vertical drawer space
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters
 - one steering caster and two braking casters

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Engineered &
Manufactured
in the U.S.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.5 cm)

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.