



# M-series

## MDS3024K14

### M-SERIES SHORT MEDICAL CART with Four Drawers and a Standard Key Lock

#### Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- 24" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.5 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



### Short Height, Standard Width, Four Drawer Medical Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered &  
Manufactured  
in the U.S.

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300