



M-series MDS3030B06+MD30-EMG3

M-SERIES TALL EMERGENCY CART with MD30-EMG3 Emergency Package

MD30-EMG3 Accessories Include:



Oxygen tank holder (O2HLDR)



Large Utility Basket (BASKETDM)



(2) Utility hooks (UHOOKS2)



Engineered & Manufactured in the U.S.



Toprail, matching cart color (MD30-TOPRAIL)



Cardiac board and brackets (MD-CARDBRD)



Chart Holder (CHARTHLDRDM)



(2 Sets) Standard 6" Drawer Dividers (MD30-DIV6-S)



(1 Set) Standard 3" Drawer Dividers (MD30-DIV3-S)



(1 Set) Standard 3" Drawer Tray Dividers (MD30-TRAYDIV3-S)



IV pole (MD-IV2)



Defibrillator Shelf (MD-DEFIB)



Locking Sharps Container (5qt) & Glove Box Holder (LKSHRPHOLD5Q-GLV)



Suction shelf (SUCTIONSHLF)

Color shown is Red. More than 20 additional colors are available.



Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE