



M-series MDS3030B06+MD30-EMG

M-SERIES TALL EMERGENCY CART with MD30-EMG Emergency Package

Tall, Standard Width, Six Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD30-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



IV pole (MD-IV2)



Cardiac board and brackets (MD-CARDBRD)

Color shown is Red. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cart is delivered assembled
- Removable plastic top
- Latex free
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- Heat sealed powder coat steel break-away lock bar
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish
- Fully customizable drawer configurations



HARLOFF

CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064
medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.