



M-series

MDS3030EKC06

M-SERIES STANDARD WIDTH CART with Six Drawers and an Electronic Lock w/Keypad

Features:

- Sturdy, uni-body steel construction
- Electronic locking cart with keypad access and Key Lock override
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- May add, edit, or delete user and user groups
- Battery operated with six D-cell Batteries; 1 year minimum life at 25-50 uses per day
- Lock has an LED low battery indicator
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- 30" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Latex free cart is delivered assembled
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



**Color shown is Navy.
More than 20 additional
colors are available.**



**Engineered &
Manufactured
in the U.S.**

Tall, Standard Width, Six Drawer Medical Cart with an Electronic Lock with Keypad

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/m-series-medical-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**