




M-series MDS3030K02-21DR

M-SERIES STANDARD WIDTH CART with Two Drawers, Storage Area and Key Lock

Includes a 21" double door storage compartment (keyed separately) and adjustable interior shelf with six possible height locations

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- 30" of vertical drawer space
- Cart is delivered assembled
- Drawer Configuration: One 3" drawer (7.6 cm), one 6" drawer (15.2 cm) and one 21" double door storage compartment (53.3 cm)
- Removable plastic top
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Fully customizable drawer configurations
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



Color shown is Navy.
More than 20 additional
colors are available.



Engineered &
Manufactured
in the U.S.

Tall, Standard Width, Two Drawer Cart and Storage Compartment with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at
<https://www.harloff.com/product-category/treatment-and-procedure-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300