



MDS3030K16+MD30-AIRWAYPKG

M-SERIES TALL TREATMENT CART with MD30-AIRWAYPKG Difficult Airway Package

MD30-AIRWAYPKG Accessories Include:



Oxygen tank holder (O2HLDR)



(2) Utility hooks (UHOOKS2)



(1) 16" Side Rail (SIDERAIL16)



Adjustable Overhead Shelf (MD30-ADJSHLF2)



5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)



4-Scope Capacity Side Scope Cabinet (680350)

Color shown is Navy. More than 20 additional colors are available.

Engineered & Manufactured in the U.S.



Features:

- Sturdy, uni-body steel construction
- Choice of lock options: www.harloff.com/locking-options/
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Aluminum side rail/push handle mounted on one side
- Industry best 5" quiet rolling casters - one steering caster and one braking casters
- Removable plastic top
- Fully customizable drawer configurations
- 30" of vertical drawer space
- Cart is latex free and delivered assembled

M-series

Drawer Configuration:

- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) & one 12" drawer (30.5 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com