




M-series

MDS3030K16+MD30-ANS3

M-SERIES TALL ANESTHESIA CART with MD30-ANS3 Deluxe Anesthesia Package

**MD30-ANS3
Accessories Include:**



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



Plastic Waste
Container
(MD-3GALWASTEDM)



Sharps Container
w/Wire Holder
(MD-SHARPS-RC)



Catheter Box w/
Rail Clip
(CATHBOX18RC)

Drawer Configuration:

- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) & one 12" drawer (30.5 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



3-Tier Back Rail
(MD30-RBR3)



(3) 12" Side Rails
(SIDERAIL12)



Large Utility Basket
(BASKET3RC)



Engineered &
Manufactured
in the U.S.



Tape and Label
Dispenser w/Rail
Clips (LABELDISP)



(1) Basic 6"
Drawer Divider
(MD30-DIV6-B)



(2) Standard 3"Drawer
Divider with Tray
(MD30-TRAYDIV3-S)

Color shown is
Navy. More than
20 additional
colors are
available.



Features:

- Sturdy, uni-body steel construction
- Choice of lock options: www.harloff.com/locking-options/
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Removable plastic top
- Fully customizable drawer configurations
- 30" of vertical drawer space
- Cart is latex free and delivered assembled



View similar products online at
<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com