



M-series

MDSS3027B14
in Stainless Steel

M-SERIES EMERGENCY CART

Stainless Steel w/Four Drawers, Breakaway Lock

*Ideal for cleaning and sterilization, this cart includes the same great M-Series features but with a stainless steel body.**

Features:

- Sturdy, uni-body stainless steel exterior construction/finish with heat sealed powder coat steel drawer pans
- Choice of lock options
- Stainless steel breakaway lock bar
- Pack of 100 breakaway seals
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- 27" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), one 9" drawer (22.9 cm) and one 12" drawer (30.4 cm)
- Latex free cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



* Not designed for use in a cart washer

Shown in stainless steel. More than 20 additional colors are available in heat sealed powder coat steel.

Stainless Steel Medium Height, Standard Width, Four Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/acute-care-carts/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300