




M-series

MDSS3027K16
in Stainless Steel

M-SERIES TALL MEDICAL CART

Stainless Steel w/Six Drawers and a Key Lock

*Ideal for cleaning and sterilization, this cart includes the same great M-Series features but with a stainless steel body.**

Features:

- Sturdy, uni-body stainless steel exterior construction/finish with heat sealed powder coat steel drawer pans
- Key lock with two keys
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- 27" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm) and three 6" drawers (15.2 cm)
- Cart is delivered assembled
- Fully curved protective bumpers
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



* Not designed for use in a cart washer

Shown in stainless steel.
More than 20 additional colors are available in heat sealed powder coat steel.

Stainless Steel Tall, Standard Width, Six Drawer Medical Cart with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300