



MPA2421E14

ALUMINUM

A-SERIES MEDIUM WIDTH CART with Four Drawers and a Standard Key Lock

Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top
- Ergonomic side push handles can be mounted on the left or right
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters
 - 1 steering caster and 2 braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 21" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown with a Pink body and drawers. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

All dimensions are available in the "MA-Dims-24wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

X-Short Height, Medium Width, Four Drawer Medical Cart with a Standard Key Lock

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at

<https://www.harloff.com/product-category/a-series-medical-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**