



Wood Style

MPA3021KC5-AN

ALUMINUM

## A-SERIES X-SHORT WOOD STYLE CART

with Five Drawers and a Standard Key Lock

### Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.

#### Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Curved drawer pulls, durable and elegant
- Key lock with two keys
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish on body, inside and out
- Industry only fully curved protective bumper guards your cart and walls from damage
- Removable plastic top, easily replaceable
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Ergonomic side push handles can be mounted on the left or right
- Pre-threaded mounting holes for accessories found at [harloff.com/product-category/accessories/](http://harloff.com/product-category/accessories/)
- 21" of vertical drawer space with fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters
  - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled



Body shown in Light Gray. More than 20 additional colors are available. Drawers shown in Wood Vinyl American Natural.

Three standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.



Engineered & Manufactured in the U.S.

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of A-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff aluminum cart.

#### Drawer Configuration:

- Three 3" drawers (7.6 cm) and two 6" drawers (15.2 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300