



MPA3030E16 with
MD30-SSFLTOP accessory


DELIVERY CART

with Stainless Steel Top + Two Drop Shelves

The cart top is designed to create a large workspace ideal for a baby scales, pediatric supplies and other equipment

Features:

- Curved drawer pulls, durable and elegant
- Stainless steel top with stainless steel drop shelves on both sides
- The large top workspace (58.5" in length when shelves are fully extended) is ideal for holding a baby scale and other large equipment
- Basic pushbutton e-lock with key lock override
- Cabinet may be programmed in the field
- Factory set user and supervisor code
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling casters - 1 steering caster and 2 braking casters
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Latex free cart is delivered fully assembled
- 30" of vertical drawer space
- Drawer Configuration: Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown with a Light Blue body and drawers. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Tall, Standard Width Delivery Cart with Stainless Steel Top and Side Flyout Shelves, Keyless Pushbutton E-Lock

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300