



Wood Style

MPA3030EC6-CM

ALUMINUM

A-SERIES TALL WOOD STYLE CART

with Six Drawers and Electronic Pushbutton Lock

Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.

Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Fully welded uni-body aluminum construction, highest quality in the industry
- Basic electronic pushbutton lock with key lock override and factory set user/supervisor codes
- Cabinet may be programmed in the field
- Curved drawer pulls, durable and elegant
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top, easily replaceable
- Pre-threaded mounting holes for accessories found at harloff.com/product-category/accessories/
- Ergonomic side push handles can be mounted on the left or right
- Fully customizable drawer configurations at no extra charge
- Cart is delivered assembled
- 30" of vertical drawer space
- Industry only fully curved protective bumper guards your cart and walls from damage
- Heat-sealed durable scratch resistant powder coat finish on body, inside and out
- Industry best 5" quiet rolling casters - one steering caster and two braking casters



Three standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.



Engineered & Manufactured in the U.S.

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of A-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff aluminum cart.

Drawer Configuration:

- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

HARLOFF
CLINICAL STORAGE BY CAPSA HEALTHCARE

View similar products online at
<https://www.harloff.com/product-category/a-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300