

# MOBILE VACCINE FREEZER CART



## COVID-19 VACCINATION FREEZER CART

### VF4300K-AC

#### Features:

- Allows for extended storage for vaccines prior to dispensing
- Freezer keeps vaccines down to  $-20^{\circ}\text{C}$
- Lab and medical-grade quality 1.8 cubic feet of freezer capacity
- Interior freezer dims are 18.5" H x 13.5" W x 13" D
- User-programmable, audible and visual high/low alarms
- Hospital grade power strip and 15 foot long cord included
- Secured drawer for storage of syringes masks, gauze and other vaccination supplies
- Key lock on drawers and freezer
- Antimicrobial plastic top, removable
- Fully welded heat sealed powder coat body and drawers with 12 year warranty

#### Colors:

Shown in Hammertone Blue, but available in more than twenty color options. View them online: [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)

Shown with Freezer,  
also available in  
Refrigerator model  
VR4300K-AC



## Mobile Vaccine Freezer Cart for the Application of COVID-19 Vaccinations, Three Drawer, Key Lock, No Accessories

Engineered & Manufactured in the U.S.



#### Warranty:

12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 51.5" H x 37.5" W x 26" D (130.8 cm H x 95.3cm W x 66cm D)
- Cart net weight: 151 lbs (68.5 kg)

#### Drawer Configuration:

- (3) 3" (7.6 cm) and (1) Compartment with pre-installed freezer

#### Optional Accessories:



Glove Box Holder  
Triple Direct Mount



Chart/document  
holders



Sharps container &  
mounting bracket



Waste container &  
mounting bracket

**HARLOFF**  
INNOVATING STORAGE FOR CARE

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300