

Mini Line Emergency Cart/ Model: 3143E

Three Drawer with Basic Electronic Pushbutton Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray,
Mauve, Red, Teal, Yellow, Forest Green, Hammertone
Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 24" (61 cm) of vertical drawer space
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 3" (7.6 cm) single wheel stem caster set, full swivel, two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) – 6" (15.2 cm)
- (2) – 9" (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5 cm)

Overall Cart Dimensions:

- 34.5" H x 19.5" W x 19.75"D
(87.6 cm H x 49.5 cm W x 50.2 cm D)
- Writing surface: 33.5" H x (81.2 cm)
- Cart net weight: 71 lbs (32.3 kg)
- Approximate shipping weight: 99 lbs (45 kg)
- Approximate shipping dims: 49" H x 24" W x 24" D
(124.5 cm H x 60.9 cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts