

## Mini Line Infection Control Cart/ Model: 3144E

Four Drawer (Short Cabinet) with Basic Electronic Pushbutton Lock, Standard Package



### Seventeen Available Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Construction:

- Painted steel cabinet and drawers

### Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 24" (61 cm) of vertical drawer space
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 3" (7.6 cm) single wheel stem caster set, full swivel, two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

### Drawer Configuration:

- (1) – 3" (7.6 cm)
- (2) – 6" (15.2 cm)
- (1) – 9" (22.9 cm)

### Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5 cm)

### Overall Cart Dimensions:

- 34.5" H x 19.5" W x 19.75"D (87.6 cm H x 49.5 cm W x 50.2 cm D)
- Writing surface: 33.5" H x (81.2 cm)
- Cart net weight: 74 lbs (33.6 kg)
- Approximate shipping weight: 101 lbs (45.9 kg)
- Approximate shipping dims: 40" H x 24" W x 24" D (101.6 cm H x 60.9 cm W x 60.9 cm D)

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts