

Mini Line Anesthesia Cart/ Model: 3156E-ANS

Six Drawer with Basic Electronic Pushbutton Lock, Specialty Accessory Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 30" (76.2 cm) of vertical drawer space
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 3" (7.6 cm) single wheel stem caster set, full swivel – two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-tier raised back rail system (MINI-RBR)
- (2) Aluminum side rails/push handles (68161)
- 6-bin tilt bin organizer (683805)

Drawer Configuration:

- (3) – 3" (7.6 cm)
- (2) – 6" (15.2 cm)
- (1) – 9" (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5 cm)

Overall Cart Dimensions:

- 56" H x 19.25" W x 19.75" D (142.2 cm H x 48.9 cm W x 50.2 cm D)
- Writing surface: 39.75" H x (101 cm)
- Cart net weight: 94 lbs (42.7 kg)
- Approximate shipping weight: 135 lbs (61.4 kg)
- Approximate shipping dims: 48" H x 24" W x 24" D (121.9 cm H x 60.9 cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts