



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini Line Anesthesia Cart with Extended Stabilizing Bumper and 5" Wheels, Specialty Accessory Package, with Key Lock #3245K-ANS

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints



Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-tier raised back rail system (MINI-RBR)
- (2) Aluminum side rails/push handles (68161)
- 6-bin tilt bin organizer (683805)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Assembled Cart Dimensions:

- 52.5" H x 19.5" W x 25" D (133.4 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 103 pounds (46.7 kg)
- Approximate shipping weight: 123 pounds (55.81 kg)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Sixteen available color choices:

***Antimicrobial Finish Available**

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White