



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini Line Phlebotomy Cart Package with Breakaway Lock #3245BPB



Features

- Stabilizing pontoon bumpers
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 24 vertical inches (60.96 cm) drawer space
- Breakaway Lock with pack of 100 numbered seals
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty

Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (GLOVEBOX1DM)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (TILTBIN6-RC)

Other Mini Line Phlebotomy Cart Options:

- Locking Bemis® Sharps Holder with 5 quart Container (WASTE3GALDM)
- Upgrade to 3 Box Glove Holder (GLOVEBOX3RC)
- Chart Holder (CHARTHLDRDM)

Specifications

Construction:

- Painted steel cabinet and drawers

Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1cm x 42.5cm)

Assembled Cart Dimensions

- 52.3" H x 27.5" W x 21.5" D
(133 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 33.8"
- Cart net weight: 108 lbs (49.1 kg)
- Approximate shipping weight: 133 pounds (28.2 kg)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>