



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Four Drawer Mini Line Phlebotomy Cart Package with Key Lock #3254KPB



Features

- Stabilizing Pontoon Bumpers
- 5-inch (12.7 cm) premium plate casters
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key-locking system with 2 keys
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in over 20 colors

Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (TILTBIN6RC)

Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (WASTE3GALDM)
- Upgrade to 3 Box Glove Holder (GLOVEBOX3RC)
- Chart Holder (CHARTHLDRDM)

Specifications

Construction:

- Painted steel cabinet and drawers

Internal Drawer Dimensions:

- 15" W x 17" D (38.1 cm W x 43.2 cm D)

Assembled Cart Dimensions

- 58.5" H x 27.5" W x 21.5" D
(148.6 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface – 40"
- Cart net weight: 128 lbs (58.2 kg)
- Approximate shipping weight: 153 pounds (69.6 kg)

Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 9-inch (22.9 cm)

Over 20 different color options available which can all be viewed here: <https://www.harloff.com/finish-options/>