

Mini Line Phlebotomy Carts/ Model: 3256LPB

Six Drawer with Turn Latch, Specialty Accessory Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- 5-inch ball bearing casters, full swivel, 2 with brake
- Cabinet has 30 vertical inches (76.2 cm) drawer space
- Turn latch to secure drawers while moving cart
- Stabilizing pontoon bumpers
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended push handle (MINI-PBHANDKIT)
- Single glove box holder (AL2362)
- Aluminum waste container with lid (AL680436)
- Two 6-bin tilt bin organizers (683805)

Other Mini Line Phlebotomy Cart Options

- Locking Bemis® sharps holder with 5 quart container (684810)
- Chart holder (680417)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (38.1 cm W x 42.5 cm D)

Overall Cart Dimensions:

- 58.5" H x 27.5" W x 21.5" D (w/ accessories in place) (148.6 cm H x 70 cm W x 54.6 cm D)
- Writing surface: 42" H (106.7 cm H)
- Cart net weight: 113 lbs (51.3 kg)
- Approximate shipping weight: 153 lbs (69.6 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (131.1 cm H x 73.7 cm W x 60.9 cm D)