

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Six Drawer Mini24 Emergency Cart (Tall Cabinet) with 3" Casters, **Breakaway Lock, Accessory Package** #4156B-EMG

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking





Additional Colors Available. • Oxygen tank holder (AL680403)

Assembled Cart Dimensions:

- 66.25" H x 30.5" W x 20.5" D (168.3 cm H x 77.5 cm W x 52 cm D) height includes IV Pole - extends to 84.75" / 215.3 cm H
- Cart net weight: 114 pounds (51.71 kg)

Shipping Dimensions:

- 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 153 pounds (69.4 kg)

Specifications

Construction:

· Painted steel cabinet and drawers

Features:

- · Space-saver design allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Breakaway lock system with package of 100 plastic seals
- 3" (7.6 cm) single-wheel stem caster set, full swivel - two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish available in seventeen colors

Accessory Package Includes:

- 2-prong extendable IV pole (680401)
- Cardiac board and brackets (680406)
- Aluminum push handle / side rail (68161)

Internal Drawer Dimensions:

• 18.5" W x 16.75" D (47 cm x 42.5cm)

Drawer Configuration:

- (3) 3-inch (7.6 cm)
- (2) 6-inch (15.2 cm)
- (1) 9-inch (22.9 cm)