



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini24 Anesthesia Cart with Bumper and 5" Wheels, with Basic Electronic Lock #4245E

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices

### Specifications



Color Shown is Sand.  
Additional Colors Available.

#### Shipping Dimensions:

- 49" H x 33" W x 26" D  
(124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 142 pounds (64.4 kg)

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) drawer space
- Basic electronic pushbutton lock with keylock override
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5" (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

#### Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

#### Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Assembled Cart Dimensions:

- 36.5" H x 29.25" W x 22" D  
(92.5 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 103 pounds (46.8 kg)