



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Six Drawer Mini24 Emergency Cart with Bumper and 5” Wheels,
 Break-away Lock
 #4256B**

Harloff’s Mini24 Line offers a 24” wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.



Color Shown is Red
 Additional Colors Available

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 30 vertical inches (76.2 cm) of drawer space
- Break-away locking system
- Industry only fully curved protective bumper
- 5” (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5” W x 16.75” D (47cm x 42.5cm)

Assembled cart dimensions:

- 42.75” H x 29.25” W x 22” D (108.4 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 119 pounds (54.1kg)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Shipping Dimensions:

- 57” H x 33” W x 26” D (144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 158 pounds (71.7 kg)