

650 Ford Street Colorado Springs, CO 80915 Toll: (800) 433-4064 * Ph: (719) 637-0300

Fax: (719) 597-8273 * www.harloff.com

Six Drawer Mini24 Anesthesia Cart with Bumper and 5" Wheels, Specialty Accessory Package, with Key Lock #4256K-ANS

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.



Color Shown is Beige.
Additional Colors Available.

Assembled Cart Dimensions:

 Cart dims: 57.75" H x 29.25" W x 22" D (147 cm H x 74.1 cm W x 56 cm D)

• Cart Net Weight: 134 pounds (60.8 kg)

Specifications

Construction:

· Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5" (12.7 cm) Premium caster set, full swivel
 two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- · Ball-bearing full extension drawer slides
- · Antimicrobial plastic top, removable
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish with seventeen available colors

Accessory Package Includes:

- 3-tier raised back rail system
- (2) Aluminum side rails/push handles (68161)
- 6-bin tilt bin organizer (683805)

Internal Drawer Dimensions:

• 18.5" W x 16.75" D (47 cm x 42.5cm)

Drawer Configuration:

- (3) 3-inch (7.6 cm)
- (2) 6-inch (15.2 cm)
- (1) 9-inch (22.9 cm)