

Stainless Steel Cast Cart/ Model: 6025E-TC

Eight Drawer with Basic Electronic Push Button Lock, Deluxe Package



Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Stainless steel cabinet and drawers
- Painted aluminum top compartment

Features:

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with vinyl mat
- Swivel tourniquet/ pump pole
- 24" stainless steel ruler etched on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Basic electronic pushbutton lock
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Key lock override
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free

Other Cast Cart Options:

- All cast carts are also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks. Call Harloff for details.

Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall Cart Dimensions:

- 43" H x 38" W x 21" D
(109.2 cm H x 96.5 cm W x 53.3 cm D)
w/ top compartment (TC) 66" H x 44" W x 21" D
(167.6 cm H x 111.7 cm W x 53.3 cm D)
- Cart net weight: 226 lbs (102.7 kg)
- Approximate shipping weight: 271 lbs (123.2 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D
(144.8 cm H x 119.4 cm W x 68.6 cm D)