



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, 2 Drawer with Lower Storage Compartment, Treatment Cart– Specialty Package #6201

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Tall cabinet
- 21” Storage compartment, locking doors and adjustable shelf
- Key lock with 4 keys
- Painted steel top rail
- Side mounted drop shelf
- 2 ea., 3” adjustable exchange trays
- 5” (12.7cm) ball-bearing premium casters, all full swivel – one with brake, one tracking
- Stabilizing pontoon bumpers
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- Storage Compartment: 21-inch (53.3 cm) high opening, with adjustable internal shelf

Internal Drawer Dimensions:

23” W x 17” D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 44.5” H x 32” W x 22” D
(113 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 114 pounds (51.7 kg)
- Approximate shipping weight: 153 pounds (69.4 kg)
- Approximate shipping dims: 53” H x 33” W x 26” D
(134.6 cm H x 83.8 cm W x 66 cm D)



Color Shown is Beige.
Additional Colors Available