

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Classic Line, Four Drawer, Emergency Cart, Standard Package #6300



Color Shown is Red. Additonal Colors Available

Specifications

Construction:

• Sturdy, uni-body steel construction

Features:

- · Painted steel break-away lock bar
- Pack of 100 break-away seals
- 24" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Stabilizing pontoon bumpers
- Five-inch (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Painted steel top rail
- Removable plastic top
- · Ball-bearing full extension drawer slides
- Cart is delivered assembled
- Twelve Year Warranty on manufactured parts
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (2) 3" (7.6 cm)
- (1) 6" (15.2 cm)
- (1) 12" (30.5 cm)

Internal Drawer Dimensions:

• 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 39.5" H x 32" W x 22" D
 - (100.3 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 104 lbs (47.2 kg)
- Approximate shipping weight: 145 lbs (65.8 kg)
- Approximate shipping dims: 47" H x 33" W x 26" D (119.4cm H x 83.8 cm W x 66.0 cm D)